

SurePath™ A Pap test

that makes sense for her.


1. Sample Collection

- ▶ Only SurePath™ can ensure 100% of the collected sample makes it to the lab every time. Just collect the sample and drop the detachable head(s) in the vial, seal and send.
- ▶ SurePath™ is the only Pap test that uses this sample collection technique.
- ▶ No swish, no swirl, no loss of diagnostically relevant cells during sample transfer¹.




2. At The Laboratory

- ▶ SurePath™ uses a proprietary Cell Enrichment™ process based on centrifugation, not filtration.
- ▶ SurePath™ Cell Enrichment™ separates and largely eliminates obscuring cellular material and debris such as mucous, inflammatory cells and blood, providing your Pathologist with better visualization of diagnostically relevant cells.



3. Results

- ▶ SurePath™ is the only Pap test available today with **a 64.4% increase in HSIL+ detection²**.
- ▶ SurePath™ significantly reduces Unsatisfactory cases². A reduction in Unsatisfactory cases help reduce unnecessary and inconvenient repeat testing.



Mary is counting on you
to choose SurePath™ for her.

SUREPATH™

liquid-based Pap test

Now You're Sure.™